

## 853-021924-003 Lam Research Paddle Lifter

### Features

- New Bellows Assembly
- New Silicon Bronze Linear Adjust Screw
- New Stepper Motor
- New Linear & Thrust Bearings
- New Thrust Washers
- New D-Sub or Viking Connectors
- Fully Disassembled
- All Parts Cleaned and Re-anodized as needed
- Helium Leak Tested
- Full Functionality Tested on an actual Lam machine



### [View Product](#)

The Lam Research Paddle Lifter cycles the heated paddle up & down with precision.

**Applications Include:** All Lam Research Standalone 4520i & 9600 Platforms.

SemiGroup offers a complete array of options and expertise when supplying Plasma Etch Equipment to our clients. We also provide Total Refurbishment Services for your existing Plasma Etch Equipment. We will refurbish all assemblies, clean and re-range all your Mass Flow Controllers, recondition all chambers and load locks, replace all o-rings and seals to return the equipment to original OEM Specifications. Your machine will go through a rigorous Plasma "On" Wafer Marathon to ensure complete functionality before we ship your tool. Any required peripheral equipment will now also be fully checked out (Pumps, Chillers, RF Carts, Cabling etc) before being shipped.

Finally, Semigroup will Startup and Install your newly refurbished etcher completing our "Turn-Key" Premium Client Service. We also have our own line of consumable parts we will supply with your newly refurbished tool.